

## GST-Tag

---

**Tested applications:**WB IF IP

**Recommended Dilution:**WB 1:2000 - 1:5000 IF 1:50 - 1:200 IP 1:20 - 1:50

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of GST tag

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**AE001

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Isotype:**IgG

**Purity:**Affinity purification

For research use only.