

## Anti-mCherry-tag Mouse mAb

---

**Tested applications:**WB

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant Protein

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**AE002

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**

For research use only.