

## TAP-tag

---

**Tested applications:**WB IP

**Recommended Dilution:**WB 1:1000 - 1:5000 IP 1:100 - 1:500

**Calculated MW:**48kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of TAP epitope

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**AE015

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Isotype:**IgG

**Purity:**Affinity purification

For research use only.