

V5-Tag

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of V5 epitope

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:AE019

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Isotype:IgG

Purity:Affinity purification

For research use only.

Background:

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids, although it has also been used with a shorter 9 amino acid sequence.

To place an order, please [Click HERE](#).