

## Soybean P34 GST

**Description:** The E.Coli derived GST Tag recombinant protein, 36kDa, contains Soybean P34 Protein epitopes 214-261, 351-379 amino acids.

Catalog #:ALPS-012

**Purity:** Protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

For research use only.

**Formulation:**

50mM Tris-HCl, pH 8.0, 60mM NaCl, 10mM glutathione and 50% glycerol.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chem.

**Introduction:**

The P34 protein is the main allergen for soybean sensitive humans. Soybean protein P34, a thiol protease belonging to the papain family, is a monomeric allergen having an N-terminal amino acid sequence and amino acid composition identical to that of the seed 34kDa protein. It is an insoluble glycoprotein having a pI of 4.5 and a calculated mass of 28.643 Dalton, representing 23% of total soybean protein. Upon glycosylation, the mass will be somewhat larger, resulting in a ~32kDa band in non-reduced SDS PAGE gels. It exhibits no enzymatic function due to an absence of the catalytic cysteine. P34 is stored in storage vacuoles of soybean cotyledons.

**Storage:**

Soybean P34 His although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

**To place an order, please [Click HERE](#).**