

Phospho-RELA-T505

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 65kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding T505 of human RELA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bkops

Synonym:

p65; NFKB3;

Catalog #: AP0190

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5970

Isotype: IgG

Swiss Prot: Q04206

Purity: Affinity purification

For research use only.

Background:

NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).