

## Phospho-ICAM1-Y512

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 89.2kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding Y512 of human ICAM1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

BB2; CD54; P3.58;

**Catalog #:** AP0365

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 3383

**Isotype:** IgG

**Swiss Prot:** P05362

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.

*To place an order, please [Click HERE](#).*