

Phospho-JUNB-S259

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 43kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S259 of human JUNB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

AP-1;

To place an order, please [Click HERE](#).

Catalog #: AP0375

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3726

Isotype: IgG

Swiss Prot: P17275

Purity: Affinity purification

For research use only.