

## Phospho-JUND-S255

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:1000 IHC 1:50 - 1:100

**Calculated MW:**38kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding S255 of human JUND

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

k

**Synonym:**

AP-1

**Catalog #:**AP0377

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3727

**Isotype:**IgG

**Swiss Prot:**P17535

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. It has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternate translation initiation site usage results in the production of different isoforms.

*To place an order, please [Click HERE](#).*