

## Phospho-Map2k5-S311/T315

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 46kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding S311/T315 of human Map2k5

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

Mek5;

**Background:**

kinase involved in signal transduction in the cell.

**To place an order, please [Click HERE](#).**

**Catalog #:** AP0392

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5607

**Isotype:** IgG

**Swiss Prot:** Q62862

**Purity:** Affinity purification

For research use only.