

## Phospho-PDPK1-S241

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:100 - 1:200

**Calculated MW:** 63kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding S241 of human PDPK1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PDK1; PDPK2; PRO0461;

*To place an order, please [Click HERE](#).*

**Catalog #:** AP0426

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5170

**Isotype:** IgG

**Swiss Prot:** O15530

**Purity:** Affinity purification

For research use only.