

Phospho-Erk1/2-T202/Y204

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 42,44KD

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding T202/Y204 of human Erk1/2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ERK; p38; p40; p41; ERK2; ERT1; MAPK2; PRKM1; PRKM2; P42MAPK; p41mapk;

Catalog #: AP0472

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5594/5595

Isotype: IgG

Swiss Prot: P27361

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the MAP kinase family. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. The activation of this kinase requires its phosphorylation by upstream kinases. Upon activation, this kinase translocates to the nucleus of the stimulated cells, where it phosphorylates nuclear targets. Two alternatively spliced transcript variants encoding the same protein, but differing in the UTRs, have been reported for this gene.

To place an order, please [Click HERE](#).