

## Phospho-KSR1-S392

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**102kDa

**Observed MW:**Refer to figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding S392 of human KSR1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

KSR; RSU2;

*To place an order, please [Click HERE](#).*

**Catalog #:**AP0487

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**8844

**Isotype:**IgG

**Swiss Prot:**Q81VT5

**Purity:**Affinity purification

For research use only.