

## Leptin Receptor Human

**Description:**Leptin Binding Domain Human Recombinant also called Leptin soluble receptor produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 208 amino acids and having a molecular mass of 24.5 kDa. Leptin Receptor consists of the cytokine binding domain of leptin receptor amino acids 428-635 of human leptin receptor.The Leptin Binding Domain is purified by proprietary chromatographic techniques.

**Catalog #:**CYPS-515

For research use only.

**Synonyms:**OB Protein, Obesity Protein, OBS, Obesity factor, Leptin Receptor.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile Filtered White lyophilized (freeze-dried) powder.

**Amino Acid Sequence:**The sequence of the first four N-terminal amino acids was determined and was found to be Ala-Thr-Pro-ValBiological ActivityBiological Activity is evidenced by high affinity binding of mammalian leptins at 1:1 molar ratio.

**Purity:**Greater than 98.0% as determined by:(a) Analysis by SEC-HPLC.(b) Analysis by SDS-PAGE.

**Formulation:**

The protein was lyophilized from a concentrated (1mg/ml) solution with 0.0045mM NaHCO<sub>3</sub>.

**Stability:**

Lyophilized Leptin Binding Domain although stable at room temperature, should be stored desiccated below 0°C. Reconstituted Leptin is best stored refrigerated at 4°C.Please prevent freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Solubility:**

The lyophilized Leptin Binding Domain is very soluble in water and most aqueous buffers below and above the isoelectric point.

**To place an order, please [Click HERE](#).**