

Acrp30 Mouse, His

Description: The Adiponectin Mouse is created as a recombinant protein with a 21 a.a N-terminal fusion of His Tag. The Adiponectin His-Tagged Fusion Protein, produced in E. coli, is a 27.2kDa protein containing 251 amino acid residues of the Acrp30 Mouse, 18-247 amino acids.

Catalog #: CYP5-544

For research use only.

Synonyms: Acrp30, AdipoQ, GBP-28, APM-1, ACDC.

Source: Escherichia Coli.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MEDDVTTEE LAPALVPPK
GTCAGWMAGI PGHPGHNGTP GRDGRDGTG EKGEKGDAGL LGPKGETGDV
GMTGAEGPRG FPGTPGRKGE PGEAAYVYRS AFSVGLETRV TVPNVPIRFT KIFYNQNH
DGSTGKFYCN IPGLYYFSYH ITVYMKDVKV SLFKKDKAVL FTYDQYQEKV VDQASGSVLL
HLEVG DQVWL QVYGDGDHNG LY

Purity: Acrp30 Mouse purity is greater than 90% as determined by SDS-PAGE.

Formulation:

Acrp30 Mouse is a sterile filtered liquid formulation containing (1mg/ml) 20mM Tris-HCl pH-8, 1mM DTT and 10% Glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Adiponectin is an adipocyte specific secreted protein that circulates in the plasma. It is induced during adipocyte differentiation and its secretion is stimulated by insulin. Mouse adiponectin shares about 83% amino acid identity with that human. Adiponectin plays a role in various physiological processes such as energy homeostasis and obesity. Adiponectin is reduced in obese humans, and decreased level is associated with insulin resistance and hyperinsulinemia.

To place an order, please [Click HERE](#).