

## BAFF Human, His

**Description:**BLyS Human Recombinant fused to His tag at N-terminus produced in E.Coli is a single, non-glycosylated polypeptide chain containing 190 amino acids and having a molecular mass of 21 kDa. The BAFF His tag fusion protein is purified by proprietary chromatographic techniques.

Catalog #:CYP5-552

For research use only.

**Synonyms:**BAFF, BLYS, CD257, TALL1, THANK, ZTNF4, TALL-1, TNFSF20, TNFSF13B, B-cell Activating Factor.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile Filtered colorless solution.

**Amino Acid Sequence:**MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHMAV  
QGPEETVTQDCLQLIADSET PTIQKGSYTF VPWLLSFKRG SALEEKENKI LVKETGYFFI  
YGQVLYTDKT YAMGHLIQRK KVHVFGDELS LVTLFRCIQN MPETLPNNSC YSAGIAKLEE  
GDELQLAIPR ENAQISLDGD VTFFGALKLL.

**Purity:**Greater than 95.0% as determined by(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.

**Formulation:**

Recombinant His Tag BAFF contains 20mM Tris buffer pH-8 and 5mM DTT.

**Stability:**

BLyS although stable at 10°C for 3 weeks, should be stored desiccated below -18°C. Please prevent freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

Binds to tnfrsf13b/taci and tnfrsf17/bcma. tnfsf13/april binds to the same 2 receptors, together, they form a 2 ligands -2 receptors pathway involved in the stimulation of b- and t-cell function and the regulation of humoral immunity. a third b-cell specific baff-receptor (baffr/br3) promotes the survival of mature b-cells and the b-cell response.B Lymphocyte Stimulator functions as a potent B-cell growth factor in costimulation assays. Administration of BAFF Human recombinant to mice disrupts splenic B-cell and T-cell zones and results in elevated levels of serum immunoglobulin.

**To place an order, please [Click HERE](#).**