

FGF12 Human

Description:The FGF-12 Human recombinant protein is a single, non-glycosylated polypeptide chain produced in E. coli, having a molecular weight of 22.6kDa and containing 201 amino acids (1-181). The FGF12 is fused to a 20 amino acid His tag at the N-terminus.

Catalog #:CYP5-627

For research use only.

Synonyms:FGF-12, FGF12, FGF12B, FHF1, Fibroblast growth factor 12, Fibroblast growth factor homologous factor 1, FHF-1, Myocyte-activating factor.

Source:Escherichia Coli.

Amino Acid Sequence:MSSHHHHHHH SSGLVPRGSH MESKEPQLKG IVTRLFSQQG
YFLQMHPDGT IDGTKDENS D YTLFNLIPVG LRVVAIQGVK ASLYVAMNGE YLYSSDVFT
PECKFKESVF ENYYVIYSST LYRQQESGRA WFLGLNKEGQ IMKGNRVKKT KPSSHFVPKP
IEVCMYREQS LHEIGEKGQR RKSSGTPTM NGGKVVNQDS T.

Purity:Greater than 90% as determined by SDS-PAGE.

Formulation:

The FGF-12 solution (1mg/ml) contains 20mM Tris pH-7.5, 1mM DTT, 2mM EDTA and 10% glycerol.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

FGF12 is part of the Fibroblast Growth Factor (FGF) family which has a vast mitogenic and cell survival functions, and play a role in a range of biological activities, among them are embryonic development, cell growth, morphogenesis, tissue repair, tumor growth, and invasion. FGF-12 doesn't obtain the N-terminal signal sequence present in the majority of the FGF family members, but it contains clusters of basic residues that act as a nuclear localization signal. When transfected into mammalian cells, FGF12 accumulated in the nucleus, but was not secreted. FGF12 is involved in nervous system development and function. FGF12 binds to IB2 (islet brain-2), a cellular kinase scaffold, and voltage gated sodium channels and is also involved in intracellular signaling and ion exchange.

Storage:

Store FGF12 at -20°C. Can be stored at 4°C for a limited period of time of 7 days.

To place an order, please [Click HERE](#).