

Dengue NS1 n

Description:The E.coli derived recombinant protein contains the NS1 Dengue Virus n-end Type-2 immunodominant regions. The Recombinant Dengue Virus NS1n is fused to a GST tag.

Catalog #:DEPS-010

Purity:Protein is >90% pure as determined by SDS-PAGE.

For research use only.

Formulation:

100mM Tris-Hcl pH 8.0, 60mM NaCl, and 5% Glycerol.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Each laboratory should determine an optimum working titer for use in its particular application.

Introduction:

Caused by one of four closely related virusserotypesof the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholinoantisense oligos have shown specific activity against Dengue virus.

Storage:

Dengue although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

To place an order, please [Click HERE](#).