

## Dengue Envelope-4 22kDa

---

**Description:**The E.coli derived recombinant 28 kDa protein is genetically engineered peptide which is derived from Dengue Type-4 Envelope. This region also contains a common antigen for Dengue type 1, 2 and 3.

Catalog #:DEPS-016

**Purity:**Protein is >95% pure as determined by 12% PAGE (coomassie staining).

For research use only.

**Formulation:**

Phosphate buffered saline, pH-7.4.

**Stability:**

Dengue Envelope ST4 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Applications:**

Each laboratory should determine an optimum working titer for use in its particular application.

**Introduction:**

Caused by one of four closely related virus serotypes of the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus.

**To place an order, please [Click HERE](#).**