

Dengue NS1, ST4

Description:The E.Coli derived recombinant protein contains the NS1 Dengue Virus full length Type-4 immunodominant regions.The dengue protein is fused to a 6xHis tag at C-terminus.

Catalog #:DEPS-024

Purity:Protein is >95% pure as determined by 12% PAGE (coomassie staining).

For research use only.

Formulation:

PBS and 0.02% NaN₃.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Each laboratory should determine an optimum working titer for use in its particular application.

Introduction:

Caused by one of four closely related virus serotypes of thegenus Flavivirus, family Flaviviridae, each serotype issufficiently different that there is no cross-protection andepidemics caused by multiple serotypes (hyperendemicity) canoccur. In cell culture experiments and mice Morpholinoantisense oligos have shown specific activity against Denguevirus.

Storage:

Dengue NS1, ST4 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

To place an order, please [Click HERE](#).