

PPP3CA Human

Description: PPP3CA Human Recombinant fused with a 23 amino acid His tag at N-terminus produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 534 amino acids (1-511 a.a.) and having a molecular mass of 60kDa. The PPP3CA is purified by proprietary chromatographic techniques.

Catalog #:ENPS-058

For research use only.

Synonyms: Serine/threonine-protein phosphatase 2B catalytic subunit alpha isoform, CAM-PRP catalytic subunit, Calmodulin-dependent calcineurin A subunit alpha isoform, PPP3CA, CALNA, CAN, CALN, CCN1, CNA1, PPP2B, CALNA1.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH TGSMSEPKAI DPKLSTTDRV
VKAVPFPSPH RLTAKEVFDN DGKPRVDILK AHLMKEGRLE ESVALRIITE GASILRQEKN
LLDIDAPVTV CGDIHQFFD LMKLFEVGGG PANTRYLFLG DYVDRGYFSI ECVLYLWALK
ILYPKTLFLL RGNHECRHLT EYFTFKQECK IKYSERVYDA CMDAFDCLPL AALMNQQFLC
VHGGLSPEIN TL

Purity: Greater than 85.0% as determined by SDS-PAGE.

Formulation:

The PPP3CA solution (0.25 mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 0.2M NaCl, 5mM DTT, 1mM EDTA and 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

PPP3CA (aka Calcineurin A) is a major soluble calmodulin binding protein in the brain and a Ca²⁺/calmodulin dependent serine/threonine protein phosphatase, with a relatively limited substrate specificity. PPP3CA activates the T cells of the immune system and can be blocked by drugs. PPP3CA activates NFATc (a transcription factor) by dephosphorylating it. The activated NFATc is subsequently translocated into the nucleus, where it upregulates the expression of interleukin 2.

To place an order, please [Click HERE](#).