

## PTP4A1 Human

**Description:** PTP4A1 Human Recombinant fused with a 20 amino acid His tag at N-terminus produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 190 amino acids (1-170 a.a.) and having a molecular mass of 21.6kDa. The PTP4A1 is purified by proprietary chromatographic techniques.

Catalog #:ENPS-059

For research use only.

**Synonyms:** Protein tyrosine phosphatase type IVA 1, PTP(CAAXI), Protein-tyrosine phosphatase 4a1, Protein-tyrosine phosphatase of regenerating liver 1, PRL-1, PTP4A1, PRL1, PTPCAAX1, HH72, DKFZp779M0721.

**Source:** Escherichia Coli.

**Physical Appearance:** Sterile Filtered colorless solution.

**Amino Acid Sequence:** MGSSHHHHHH SSGLVPRGSH MARMNRPAPV EVTYKNMRF  
ITHNPTNATL NKFIEELKKY GVTTIVRVCE ATYDTTLVEK EGIHVLDWPF DDGAPPSNQI  
VDDWLSLVKI KFREEPGCCI AVHCVAGLGR APVLVALALI EGGMKYEDAV QFIRQKRRGA  
FNSKQLLYLE KYRPMRLRF KDSNGHRNNC.

**Purity:** Greater than 90.0% as determined by SDS-PAGE.

### Formulation:

The PTP4A1 solution (1 mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 1mM DTT.

### Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

### Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Introduction:

Protein tyrosine phosphatase type IVA 1 (PTP4A1) is a cell signaling molecule that has regulatory roles in a variety of cellular processes. PTP4A1 is a unique nuclear PTP which is induced in regenerating liver and mitogen stimulated cells. PTP4A1 is primarily expressed in the spleen, bone marrow, thymus, lymph nodes, T lymphocytes and tonsil and is overexpressed in tumor cell lines.

**To place an order, please [Click HERE](#).**