

## TRPT1 Human

**Description:**TRPT1 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 277 amino acids (1-253) and having a molecular mass of 30.3kDa. TRPT1 is fused to a 24 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:ENPS-215

For research use only.

**Synonyms:**tRNA 2'-phosphotransferase 1, TRPT1.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile Filtered colorless solution.

**Amino Acid Sequence:**MGSSHHHHHH SSGLVPRGSH MGSHMNFSGG GRQEAAGSRG  
RRAPRPREQD RDVQLSKALS YALRHGALKL GLPMGADGFV PLGTLQLLPQ FRGFSAEDVQ  
RVVDNRKQR FALQLGDPST GLLIRANQGH SLQVPKLELM PLETPQALPP MLVHGTFWKH  
WPSILLKGLS CQGRTHIHLA PGLPGDPGII SGMRSHEIA VFIDGPLALA DGIPFFRSAN  
GVILTPGNTD GF

**Purity:**Greater than 95.0% as determined by SDS-PAGE.

### Formulation:

The TRPT1 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 1mM DTT.

### Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

### Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Introduction:

tRNA 2'-phosphotransferase 1 (TRPT1) is a member of the KptA/TPT1 family. TRPT1 catalyzes the last step of tRNA splicing, the transfer of the splice junction 2'-phosphate from ligated tRNA to NAD to produce ADP-ribose 1''-2'' cyclic phosphate. TRPT1 is broadly expressed; it is weakly or not expressed in the lung, spleen, small intestine and peripheral blood leukocytes.

**To place an order, please [Click HERE](#).**