

DLD Human

Description: DLD Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 511 amino acids (36-509 a.a.) and having a molecular mass of 54.4 kDa. The DLD is fused to a 37 amino acid His Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: ENPS-509

For research use only.

Synonyms: EC 1.8.1.4, DLD, DLDH, GCSL, PHE3, Dihydrolipoyl dehydrogenase mitochondrial, Dihydrolipoamide dehydrogenase, Glycine cleavage system L protein, LAD, E3.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear colorless solution.

Amino Acid Sequence: MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMADQ
PIDADVTVIG SGPGGYVAAI KAAQLGFKTV CIEKNETLGG TCLNVGCIPS KALLNNSHYY
HMAHGKDFAS RGIEMSEVRL NLDKMMEQKS TAVKALTGGI AHLFKQNKVV HVNGYGKITG
KNQVTATKAD GGTQVIDTKN ILIATGSEVT PFPGITDED TIVSSTGALS LKKVPEKMVV
IGAGVIGVEL GS

Purity: Greater than 95% as determined by SDS-PAGE.

Formulation:

The DLD solution contains 20mM Tris-HCl pH-8, 1mM DTT, 0.1M NaCl and 10% glycerol.

Stability:

DLD Recombinant Human although stable at 4°C for 30 days, should be stored below -20°C for periods greater than 30 days. Please avoid freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

DLD is an L protein of the mitochondrial glycine cleavage system which is also a component of the pyruvate dehydrogenase complex, the alpha-ketoglutarate dehydrogenase complex, and the branched-chain alpha-keto acid dehydrogenase complex. DLD mutations were found in patients with E3-deficient maple syrup urine disease and lipoamide dehydrogenase deficiency.

To place an order, please [Click HERE](#).