

UBE2V2 Human

Description:UBE2V2 Recombinant Human produced in E.Coli is a single, non-glycosylated polypeptide chain containing 165 amino acids (1-145 a.a.) and having a molecular mass of 18.5 kDa. The UBE2V2 is fused to 20 amino acid His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #:ENPS-557

For research use only.

Synonyms:DDVit-1, DDVIT1, EDAF-1, EDPF-1, EDPF1, MMS2, UEV-2, UEV2, Ubiquitin-conjugating enzyme E2 variant 2, Enterocyte differentiation-associated factor 1, Enterocyte differentiation-promoting factor 1, Vitamin D3-inducible protein, UBE2V2.

Source:Escherichia Coli.

Physical Appearance:Sterile filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MAVSTGVKVP RNFRLLEELE
EGQKGVGDGT VSWGLEDDED MTLTRWTGMI IGPPRTNYEN RIYSLKVECG PKYPEAPPSV
RFVTKINMNG INNSSGMVDA RSIPVLAKWQ NSYSIKVVLQ ELRRLMMSKE NMKLPQPPEG
QTYNN.

Purity:Greater than 95.0% as determined by SDS-PAGE.

Formulation:

UBE2V2 Human solution containing 20mM Tris pH-8, 0.1M NaCl & 10% glycerol.

Stability:

UBE2V2 Human although stable at 4°C for 1 week, should be stored desiccated below -18°C.
Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

UBE2V2 comprises a distinct subfamily within the E2 protein family. They have sequence similarity to other ubiquitin-conjugating enzymes however require the conserved cysteine residue that is vital for the catalytic activity of E2s. UBE2V2 shares homology with ubiquitin-conjugating enzyme E2 variant 1 and yeast MMS2 gene product. UBE2V2 participates in the differentiation of monocytes and enterocytes.

To place an order, please [Click HERE](#).