

GMPR Human

Description: GMPR Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 365 amino acids (1-345) and having a molecular mass of 39.5kDa. GMPR is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: ENPS-577

For research use only.

Synonyms: GMP reductase 1, Guanosine 5'-monophosphate oxidoreductase 1, Guanosine monophosphate reductase 1, GMPR, GMPR1.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MPRIDADLKL DFKDVLLRPK
RSSLSRAEV DLERTFTFRN SKQTYSGIPI IVANMDTVGT FEMAAMVMSQH SMFTAIHKHY
SLDDWKL FAT NHPECLQNVV VSSGSGQNDL EKMTSILEAV PQVKFICLDV ANGYSEHFVE
FVKLVRAKFP EHTIMAGNVV TGEMVEELIL SGADIIVKGV GPGSVCTTTRT KTGVGYPQLS
AVIECADSAH GL

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

The GMPR solution (1mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 40% glycerol, 0.15M NaCl and 1mM DTT.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Guanosine monophosphate reductase (GMPR) catalyzes the irreversible NADPH-dependent deamination of GMP to IMP. GMPR acts in the conversion of nucleobase, nucleoside and nucleotide derivatives of G to A nucleotides, and in upholding the intracellular balance of A and G nucleotides. In addition, the GMPR protein functions in the re-utilization of free intracellular bases and purine nucleosides.

To place an order, please [Click HERE](#).