

## HCV Mosaic-A

---

**Description:**HCV Mosaic-A protein contains a short core peptide residue 1-50a.a., 1269-1316a.a. of NS3, three epitopes from NS4 and two epitopes from NS5 and the genotype of all sequences are 1b. A specific peptide was identified from each region; the four peptides from each region were linked together and expressed in E .coli. The recombinant protein migrates at 33kDa.

**Catalog #:**HCPS-008

For research use only.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile Filtered clear solution.

**Purity:**Protein is >95% pure as determined by 12% PAGE (coomassie staining).

**Formulation:**

HCV Mosaic-A protein solution containing PBS.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Applications:**

ELISA, gold conjugation.

**Introduction:**

Hepatitis C is one type of viral hepatitis - a liver disease - caused by the hepatitis C virus (HCV) which is usually passed through contact with infected blood but can also pass through sex with an infected person and from mother to baby during childbirth.HCV core protein, NS3, NS4, and NS5 proteins of hepatitis C virus are all essential for detection of antibodies against HCV.

**Storage:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles.

**To place an order, please [Click HERE](#).**