

HSP27 Human, His

Description: Recombinant Human HSP-27 produced in E.Coli is a single, non-glycosylated polypeptide chain containing 225 amino acids (1-205 a.a.) and having a molecular mass of 24.9kDa. Hsp27 is fused to a 20 amino acid his tag at N-terminus and purified by conventional chromatography.

Catalog #: HYP5-036

For research use only.

Synonyms: HSPB1, CMT2F, HMN2B, HSP27, HSP28, HSP25, Heat shock protein beta-1, Heat shock 27 kDa protein, Stress-responsive protein 27, SRP27, HS.76067, DKFZp586P1322.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MTERRVPFSL LRGPSWDPFR
DWYPHSRLFD QAFGLPRLPE EWSQWLGGSS WPGYVRPLPP AAIESPAVAA PAYSRALSRQ
LSSGVSEIRH TADRWRVSLD VNHFAPELDT VKTKDGVVEI TGKHEERQDE HGYISRCFTR
KYTLPPGVDP TQVSSLSPE GTLTVEAPMP KLATQSNEIT IPVTFESRAQ LGGPEAAKSD
ETAAK.

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

The HSP27 protein solution contains 20mM Tris-Hcl buffer, 1mM DTT and 10% glycerol.

Stability:

HSP-27 his tag although stable at 4°C for 1 week, should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HSP27 (also known as the Estrogen-Regulated 24K protein, and hsp 28) is a member of the mammalian small heat shock protein family. Hsp27 is expressed constitutively in many tissues and its expression is increased to high levels after various types of stress including elevated temperatures, toxic metals, drugs and oxidants. Also, Hsp27 is phosphorylated in vivo on three phosphorylation sites (Ser15, Ser78 and Ser82) by protein kinases including MAPKAP kinase 2 and the stress-activated protein kinase SAPK2 (p38).

To place an order, please [Click HERE](#).