

HSPB9 Human

Description: HSPB9 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 179 amino acids (1-159 a.a.) and having a molecular mass of 19.6kDa. HSPB9 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: HYP5-058

For research use only.

Synonyms: Heat shock protein beta-9, HspB9, Cancer/testis antigen 51, CT51, HSPB9.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MQRVGNTFSN ESRVASRCPS
VGLAERNRVA TMPVRLLRDS PAAQEDNDHA RDGFQMKLDA HGFAPPEELVV QVDGQWLMVT
GQQQLDVRDP ERVSYRMSQK VHRKMLPSNL SPTAMTCCLT PSGQLWVRGQ CVALALPEAQ
TGSPRLGSL GSKASNLTR.

Purity: Greater than 85% as determined by SDS-PAGE.

Formulation:

The HSPB9 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 50% glycerol, 0.2M NaCl and 2mM DTT.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Heat Shock Protein B9 (HSPB9) is a member of the small heat shock protein (HSP20) family. HSPB9 is found in the cytoplasm as well as in the nucleus. HSPB9 is expressed specifically in the testis.

To place an order, please [Click HERE](#).