

Monkey MPIF-1 ELISA Kit

Description:Monkey Myeloid Progenitor Inhibitory Factor 1

SAMPLE TYPE:Cell culture supernatant, Urine, Serum, Plasma

CROSS REACTS WITH:Monkey

PRODUCT OVERVIEW:The ELISA (Enzyme-Linked Immunosorbent Assay) kit is an in vitro enzyme-linked immunosorbent assay for the quantitative measurement of samples in serum, plasma, cell culture supernatants and urine.

INTENDED USE:This Monkey MPIF-1 ELISA kit is intended for laboratory research use only and not for use in diagnostic or therapeutic procedures. The stop solution changes color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Monkey MPIF-1 in the sample, this Monkey MPIF-1 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of optical density versus Monkey MPIF-1 concentration. The concentration of the samples is then determined by comparing the O.D. of the samples to the standard curve.

STORAGE INSTRUCTIONS:Store at 4°C

SAFETY NOTES:

1. This kit contains materials with small quantities of sodium azide. Sodium azide reacts with lead and copper plumbing to form explosive metal azides. Upon disposal, flush drains with a large volume of water to prevent azide accumulation. Avoid ingestion and contact with eyes, skin or mucous membranes. In the case of contact, rinse the affected area with plenty of water. Observe all federal, state and local regulations for disposal.
2. All blood components and biological materials should be handled as potentially hazardous. Follow universal precautions as established by the Centers for Disease Control and Prevention and by the Occupational Safety and Health Administration when handling and disposing of infectious agents.

To place an order, please [Click HERE](#).

Catalog #:MKM0342

Sample Volume:100 ul

Assay Time:90 minutes

Principle:Competitive ELISA

Instrument:Microplate Reader

Conjugate:HRP

Detection Method:Colorimetric

COMPONENTS:

Microtiter Plate

Enzyme Conjugate

Standards

Substrate A

Substrate B

Stop Solution

Wash Solution

Balance Solution

For research use only.