

CCNB2 Human

Description:CCNB2 Human Recombinant produced in E. coli is a single polypeptide chain containing 422 amino acids (1-398) and having a molecular mass of 47.9 kDa.CCNB2 is fused to a 24 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PKPS-042

Synonyms:G2/mitotic-specific cyclin-B2, HsT17299, cyclin B2.

For research use only.

Source:E.coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MGSHMALLRR PTVSSDLENI
DTGVNSKVKS HVTIRRTVLE EIGNRVTTTRA AQVAKKAQNT KVPVQPTKTT NVNKQLKPTA
SVKPVQMEKL APKGPSPTPE DVSMKEENLC QAFSDALLCK IEDIDNEDWE NPQLCSDYVK
DIYQYLRQLE VLQSINPHFL DGRDINGRMR AILVDWLVQV HSKFRLLQET LYMCVGIMDR
FLQVQPVSRK KL

Purity:Greater than 90% as determined by SDS-PAGE.

Formulation:

The CCNB2 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.2M NaCl, 5mM DTT and 50% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

CCNB2 is a member of the cyclin family. CCNB2 is vital for regulation of the cell cycle at the G2/M (mitosis) transition. CCNB2 cooperates with the CDK1 protein kinase to create a serine/threonine kinase holoenzyme complex recognized as maturation promoting factor (MPF).

To place an order, please [Click HERE](#).