

## CK2h Human

---

**Description:** Human recombinant casein kinase 2 holo enzyme containing alpha and beta subunits which were separately expressed in E. coli as non-fusion proteins and purified using several chromatographic steps. The holo enzyme has been reconstituted in the course of the purification and is highly active suitable for labelling CK2 substrates. CK2 holoenzyme is a non-glycosylated polypeptide having a molecular mass of 140 kDa.

**Catalog #:** PKPS-218

For research use only.

**Synonyms:** Casein Kinase 2 Holoenzyme, CK2h.

**Source:** Escherichia Coli.

**Physical Appearance:** Sterile Filtered clear solution.

**Purity:** Greater than 95% as determined by SDS-PAGE.

**Formulation:**

CK2 Holoenzyme is supplied 0.5mg/1ml in 25mM Tris-HCl, 500mM NaCl, 1mM DTT, 500

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Storage:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

**To place an order, please [Click HERE](#).**