

PKC-t Human

Description:PKC-q Human Recombinant produced in Sf9 is a glycosylated, polypeptide chain containing amino acids 2-706 and having a molecular mass of 86 KD. This protein is the full-length form of the protein with an amino terminal poly His-tag.PKC-q is purified by proprietary chromatographic techniques.

Catalog #:PKPS-331

For research use only.

Synonyms:Protein kinase C theta type, EC 2.7.11.13, nPKC-theta, PRKCQ, PRKCT, MGC126514, MGC141919.

Source:Sf9 insect cells.

Physical Appearance:Sterile Filtered clear solution.

Purity:Greater than 85% as determined by SDS-PAGE.

Formulation:

PKC-q is supplied at a of 0.1mg/ml in 10mM Tris, pH 7.4, 0.1M NaCl, 20% glycerol, 1mM DTT, 0.1mM EDTA, 0.2mM PMSF and 0.03% Brij-35.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play a distinct role. The protein encoded by this gene is one of the PKC family members. It is a calcium-independent and phospholipid-dependent protein kinase. This kinase is important for T-cell activation. It is required for the activation of the transcription factors NF-kappaB and AP-1, and may link the T cell receptor (TCR) signaling complex to the activation of the transcription factors.

Storage:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

To place an order, please [Click HERE](#).