

## PI3KD946GST Human

**Description:**The PI3KD946A protein is a catalytically inactive mutant of PI3K with a D946A mutation in the ATP binding site. This recombinant catalytically inactive protein can be used as a negative control in any kind of PI3K kinase activity studies. Recombinant full length PI3KgDA mutant (126.3 kDa without Tag) carries an N-terminal GST-Tag which facilitates the proteins application in typical GST pull-down assays.

Catalog #:PKPS-338

For research use only.

**Synonyms:**Phosphoinositide 3-kinase p110, PI3KD946.

**Source:**Sf9 insect cells.

**Physical Appearance:**Sterile filtered liquid formulation.

**Purity:**Greater than 95.0% as determined by SDS Page.

**Formulation:**

0.4mg/ml solution in 25mM Tris-HCl, pH 8.0, 50mM NaCl, 0.5mM MgCl<sub>2</sub>, and 50% glycerol.

**Stability:**

PI3Kb although stable at 4°C for 1 week, should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Please prevent freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**To place an order, please [Click HERE](#).**