

VAT1 Human

Description:VAT1 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 413 amino acids (1-393 a.a.) and having a molecular mass of 44.1kDa. VAT1 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-1018

For research use only.

Synonyms:Synaptic vesicle membrane protein VAT-1 homolog, VAT1, VATI.

Source:E.coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MSDEREVAEA ATGEDASSPP
PKTEAASDPQ HPAASEGAAA AAASPPLLRC LVLTFGGYD KVKLQSRPAA PPAPGPGQLT
LRLRACGLNF ADLMARQLY DRLPPLPVTP GMEGAGVVIA VEGEVSDRKA GDRVMVLNRS
GMWQEEVTVP SVQTFLLPEA MTFEEAAALL VNYITAYMVL FDFGNLQPGH SVLVHMAAGG
VGMAAVQLCR T

Purity:Greater than 90% as determined by SDS-PAGE.

Formulation:

VAT1 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 10% glycerol and 100mM NaCl.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Synaptic vesicle membrane protein VAT-1 homolog (VAT1) is a member of the quinone oxidoreductase subfamily of zinc-containing alcohol dehydrogenase proteins. Synaptic vesicles are in charge of regulating the storage and release of neurotransmitters in the nerve terminal. VAT1 has an increased calcium ion-dependent expression in glioblastomas and on wounding, in basal keratinocytes. VAT1 is an abundant integral membrane protein of cholinergic synaptic vesicles and is believed to be involved in vesicular transport.

To place an order, please [Click HERE](#).