

LRPAP1 Human

Description:LRPAP1 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 348 amino acids (35-357) and having a molecular mass of 40.4kDa.LRPAP1 is fused to a 25 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-1036

For research use only.

Synonyms:Low density lipoprotein-related protein-associated protein 1 (alpha-2-macroglobulin receptor-associated protein 1), A2MRAP, HBP44, RAP, A2RAP, Alpha-2-MRAP.

Source:E.coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MGSHMYSREK NQPKPSPKRE
SGEEFRMEKL NQLWEKAQRL HLPPVRLAEL HADLKIQRD ELAWKKLKLD GLDEDEGEKEA
RLIRNLNVIL AKYGLDGKGD ARQVTSNSLS GTQEDGLDDP RLEKLWHKAK TSGKFSGEEL
DKLWREFLHH KEKVHEYNLV LETLSRTEEI HENVISPSDL SDIKGSVLHS RHTELKEKLR
SINQGLDRLR RV

Purity:Greater than 90% as determined by SDS-PAGE.

Formulation:

The LRPAP1 solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT and 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

LRPAP1 is a member of the alpha-2-MRAP family. LRPAP1 cooperates with LRP1/alpha-2-macroglobulin receptor and glycoprotein 330. LRPAP1, located on the cell surface, creates a complex with the alpha 2 macroglobulin receptor light and heavy chains and has part in the pathogenesis of membrane glomerular nephritis.

To place an order, please [Click HERE](#).