

## ZFAND3 Human

**Description:**ZFAND3 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 251 amino acids (1-227) and having a molecular mass of 27.7kDa.ZFAND3 is fused to a 24 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

**Catalog #:**PRPS-1037

For research use only.

**Synonyms:**Zinc finger, AN1-type domain 3, TEX27, Testis-expressed sequence 27, FLJ13222, AN1-type zinc finger protein 3.

**Source:**E.coli.

**Physical Appearance:**Sterile Filtered colorless solution.

**Amino Acid Sequence:**MGSSHHHHHH SSGLVPRGSH MGSHMGDAGS ERSKAPSLPP  
RCPGFWGSS KTMNLCSKCF ADFQKKQPDD DSAPSTSNSQ SDLFSEETTS DNNNTSITTP  
TLSPSQQLP TELNVTSPSK EECGPCTDTA HVSLITPTKR SCGTDSQSEN EASPVKRPRL  
LENTERSEET SRSKQKSRRR CFQCQTKLEL VQQLGSCRC GYVFCMLHRL PEQHDCTFDH  
MGRGREEAIM KM

**Purity:**Greater than 90% as determined by SDS-PAGE.

### Formulation:

The ZFAND3 solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 100mM NaCl, 1mM DTT and 30% glycerol.

### Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

### Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Introduction:

ZFAND3 holds a DNA-binding domain and has a large number of purposes, most of which include certain form of transcriptional activation or repression. ZFAND3 is a 251 aa protein having two AN1-type zinc fingers and two UIM (ubiquitin-interacting motif) repeats. The AN1-type zinc finger domain is conserved in animals and plants and is frequently found in proteins which hold an ubiquitin-like domain, which proposes a part in the ubiquitination pathway.

**To place an order, please [Click HERE](#).**