

SEC22B Human

Description:SEC22B Human Recombinant produced in E. coli is a single polypeptide chain containing 205 amino acids (14-194) and having a molecular mass of 23.3 kDa. SEC22B is fused to a 24 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-1219

For research use only.

Synonyms:SEC22 Vesicle Trafficking Protein Homolog B (S. Cerevisiae) Gene/Pseudogene, SEC22 Vesicle Trafficking Protein-Like 1 (S. Cerevisiae), ER-Golgi SNARE Of 24 KDa, SEC22L1, ERS-24.

Source:E.coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MGSM LPLAAS MQEDEQSGRD
LQQYQSQAKQ LFRKLNQSP TRCTLEAGAM TFHYIEQGV CYLVLCEAAF PKKLAFAYLE
DLHSEFDEQH GKKVPTVSRP YSFIEFDTFI QKTKKLYIDS RARRNLGSIN TELQDVQRIM
VANIEEVLQR GEALSALDSK ANNLSLSSKK YRQDAKYLNM RSTYA

Purity:Greater than 90% as determined by SDS-PAGE.

Formulation:

The SEC22B solution (0.5mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 2M Urea, 0.15M NaCl and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Vesicle-trafficking protein SEC22b precursor (SEC22B) is a member of the synaptobrevin family. SEC22B forms a complex with SNARE and it is believed to have a role in the ER-Golgi protein trafficking. SEC22B has a strong similarity to the mouse and Chinese hamster proteins.

To place an order, please [Click HERE](#).