

RWDD1 Human

Description: RWDD1 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 266 amino acids (1-243 a.a.) and having a molecular mass of 30.3kDa (Molecular size on SDS-PAGE will appear higher). RWDD1 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-1316

For research use only.

Synonyms: RWD domain-containing protein 1, DRG family-regulatory protein 2, RWDD1, DFRP2, CGI-24, PTD013.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGSMTDYGEE QRNELEALES
IYPDSFTVLS ENPPSFTITV TSEAGENDET VQTTLKFTYS EKYPDEAPLY EIFSQENLED
NDVSDILKLL ALQAEENLGM VMIFTLVTAV QEKLNEIVDQ IKTRREEEEKK QKEKEAEEAE
KQLFHGTPVT IENFLNWKAK FDAELLEIKK KRMKEEQAG KNKLSGKQLF ETDHNLDTSD
IQFLEDAGNN VE

Purity: Greater than 85.0% as determined by SDS-PAGE.

Formulation:

RWDD1 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.1M NaCl, 10% glycerol and 2mM DTT.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

RWD Domain Containing 1 (RWDD1) is a member of the RWDD1/GIR2 family. The RWDD1 protein protects DRG2 (developmentally regulated GTP binding protein 2), from proteolytic degradation. RWD domain refers to 3 signature motifs in proteins: RING finger, WD-repeats, and yeast DEAD (DEXD)-like motif. Current evidence proposes that RWD domain may be required for protein-protein interaction.

To place an order, please [Click HERE](#).