

BPI Human

Description:The native protein BPI is isolated from human buffy coats and has a molecular mass of 53kDa.

Catalog #:PRPS-136

Synonyms:Bactericidal/permeability-increasing protein, CAP 57.

For research use only.

Purity:Greater than 70% as determined by SDS-PAGE and capillary electrophoresis.

Formulation:

BPI (1.15mg/ml) is supplied in 25mM Sodium Acetate pH-4.6.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

BPI encodes a lipopolysaccharide binding protein. BPI is related to human neutrophil granules and has bactericidal activity on gram-negative organisms.

Storage:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

To place an order, please [Click HERE](#).