

TUBB3 Human

Description:TUBB3 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 473 amino acids (1-450 a.a.) and having a molecular mass of 52.8kDa. TUBB3 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-1395

For research use only.

Synonyms:Beta-4, CDCBM, CFEOM3A, TUBB4 , Tubulin beta-3 chain, Tubulin beta-4 chain, Tubulin beta-III, TUBB3.

Source:E.coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MGSMREIVHI QAGQCGNQIG
AKFWEVISDE HGIDPSGNYV GSDLQLERI SVYYNEASSH KYVPAILVD LEPGTMDSVR
SGAFGHLFRP DNFIFGQSGA GNNWAKGHYT EGAELVDSVL DVVRKECENC DCLQGFQLTH
SLGGGTGSGM GTLLISKVRE EYPDRIMNTF SVVSPKVSD TVVEPYNATL SIHQLVENTD
ETYCIDNEAL YD

Purity:Greater than 90% as determined by SDS-PAGE.

Formulation:

TUBB3 protein solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 10% glycerol and 1mM DTT.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Nuclear Receptor-Interacting Protein (NRIP3) is a 241 amino acid protein which is encoded by a gene that maps to human chromosome 11. Diseases associated with NRIP3 include sarcoma, and acute myeloid leukemia.

To place an order, please [Click HERE](#).