

UBL5 Human

Description:UBL5 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 93 amino acids (1-73a.a) and having a molecular mass of 10.7kDa.UBL3 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-238

For research use only.

Synonyms:Ubiquitin-like protein 5, HUB1, FLJ46917, MGC131795.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered clear solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MIEVVCNDRL GKKVRVKCNT
DDTIGDLKKL IAAQTGTRWN KIVLKKWYTI FKDHVSLGDY EIHDGMNLEL YYQ

Purity:Greater than 90.0% as determined by SDS-PAGE.

Formulation:

UBL5 protein solution (1mg/1ml) containing 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 2mM DTT.

Stability:

UBL5 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

UBL5 is believed to be a reversible modulator of protein function unlike ubiquitin which is a protein degrader. UBL5 influences a nuclear phase obligatory for increasing a reaction against the risk of mitochondrial protein misfolding.

To place an order, please [Click HERE](#).