

SCG5 Human

Description: SCG5 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 195 amino acids (27-212 a.a.) and having a molecular mass of 22 kDa. SCG5 is fused to a 8 amino acid His tag at C-terminus and purified by standard chromatography techniques.

Catalog #: PRPS-474

For research use only.

Synonyms: 7B2, SgV, P7B2, SGNE1, Secretogranin-5, Secretory granule endocrine protein I.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless liquid formulation.

Amino Acid Sequence: MYSRTPDRV SEADIQRLLH GVMEQLGIAR PRVEYPAHQA
MNLVGPQSIE GGAHEGLQHL GPFGNIPNIV AELTGDNIPK DFSEDQGYDPPNCPVGKT
ADDGCLNTPT DTAEFSREFQ LHQHLSDEPH DYPGLGKWNK KLLYEKMKGG ERRKRRSVNP
YLQGQRLDNV VAKKSVPHFS DEDKDPELEH HHHHH.

Purity: Greater than 85% as determined by SDS-PAGE.

Formulation:

FABP3 His-Tag 0.5mg/ml protein solution is supplied in 20mM Tris-HCl pH-8, 0.1mM PMSF, 1mM DTT, 0.1M NaCl and 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

SCG5 also called 7B2 protein is localized in neuroendocrine tissues. SCG5 functions as a chaperone protein for the proprotein convertase PC2 and is essential for the production of an active PC2 enzyme.

To place an order, please [Click HERE](#).