

CHGA Human, GST

Description:Chromogranin A Human Recombinant (260-454 amino acids) expressed in E.coli, shows a 47kDa band on SDS-PAGE. The Chromogranin A is fused to a 26kDa GST-tag and purified by proprietary chromatographic techniques.

Catalog #:PRPS-498

For research use only.

Synonyms:Chromogranin-A, CgA, Pituitary secretory protein I, SP-I, CHGA.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered clear solution.

Formulation:

Chromogranin A at 100

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months.Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

ELISA Inhibition Assays Western Blotting.

Introduction:

Chromogranin A (a protein of 439 amino acid which is encoded on chromosome 14) is present in neuroendocrine cells throughout the body, including the neuroendocrine cells of the large and small intestine, adrenal medulla and pancreatic islets. It is an excellent marker for carcinoid tumors, pheochromocytomas, paragangliomas, and other neuroendocrine tumors.Coexpression of chromogranin A and neuron specific Neuron Specific Enolase (NSE) is common in neuroendocrine neoplasms.

To place an order, please [Click HERE](#).