

## ANXA8 Human

**Description:**ANXA8 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 347 amino acids (1-327 a.a.) and having a total molecular mass of 39 kDa. ANXA8 is fused to a 20 amino acid His Tag at N-terminus and is purified by proprietary chromatographic techniques.

Catalog #:PRPS-555

For research use only.

**Synonyms:**ANX8, Annexin-8, Annexin VII, Vascular anticoagulant-beta, VAC-beta.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile Filtered colorless clear solution.

**Amino Acid Sequence:**MGSSHHHHHH SSGLVPRGSH MAWWKSWIEQ EGVTVKSSSH  
FNPDPDAETL YKAMKGIGTN EQAIIDVLTK RSNTQRQQIA KSFKAQFGKD LTETLKSELS  
GKFERLIVAL MYPPYRYEAK ELHDAMKGLG TKEGVIIIEIL ASRTKNQLRE IMKAYEEDYG  
SSLEEDIQAD TSGYLERILV CLLQGSRDDV SSFVDPGLAL QDAQDLYAAG EKIRGTDEMK  
FITILCTRSA TH

**Purity:**Greater than 90.0% as determined by SDS-PAGE.

### Formulation:

The ANXA8 1mg/ml protein solution contains 20mM Tris-HCL pH-8, 1mM DTT, 1mM EDTA, 0.1M NaCl and 10% Glycerol.

### Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.

### Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Introduction:

ANXA8 is part of the annexin family of evolutionarily conserved Ca<sup>2+</sup> and phospholipid binding proteins. ANXA8 is an anticoagulant protein that functions as an indirect inhibitor of the thromboplastin-specific complex, which takes part in the blood coagulation cascade. Where co-expressed in the same tissues, ANXA8 is expressed at a 100-fold lower level than Annexin A5. Though, ANXA8 is preferentially expressed in acute promyelocytic leukemia (APL) cells which suggests its involvement in hematopoietic cell differentiation.

**To place an order, please [Click HERE](#).**