

C1QTNF1 Human

Description: C1QTNF1 Human Recombinant produced in E.Coli is a single, non-glycosylated, Polypeptide chain containing 266 amino acids and having a molecular mass of 30.45 kDa. The protein contains an extra His tag at N-terminus. The C1QTNF1 amino acid sequence is identical to UniProtKB/Swiss-Prot entry Q9BXJ1 amino acids 26281. The C1QTNF1 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-661

For research use only.

Synonyms: GIP, CTRP1, ZSIG37, FLJ90694, C1QTNF1, Complement C1q Tumor Necrosis Factor-Related Protein 1, G protein-coupled receptor-interacting protein.

Source: Escherichia Coli.

Amino Acid Sequence: MKHHHHHHAS RVPHVQGEQQ EWEGTEELPS PPDHAERAE
QHEKYRPSQD QGLPASRCLR CCDPGTSMYP ATAVPQINIT ILKGEKGDRG DRGLQGKYGK
TGSAGARGHT GPKGQKGSMD APGERCKSHY AAFSVGRKKP MHSNHYYQTV IFDTEFVNLY
DHFNMFTGKF YCYVPGLYFF SLNVHTWNQK ETYLMHMKNE EEVILFAQV GDRSIMQSQS
LMLELREQDQ VW

Purity: Greater than 95% as determined by SDS PAGE.

Formulation:

Human C1QTNF1 was lyophilized from 50mM Acetate Buffer pH-4.

Stability:

Store lyophilized protein at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles. Reconstituted protein can be stored at 4°C for a limited period of time; it does not show any change after two weeks at 4°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

Add 0.1M Acetate buffer pH4 to prepare a working stock solution of approximately 0.5mg/mL and let the lyophilized pellet dissolve completely. For conversion into higher pH value, we recommend intensive dilution by relevant buffer to a concentration of 10g/ml. In higher concentrations the solubility of this antigen is limited. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

Introduction:

C1QTNF1 is a novel adipokine, providing a significant framework to further address the physiological functions and mechanisms of the action of this family of secreted glycoproteins in normal and disease states. C1QTNF1 increases the production of aldosterone. C1QTNF1 is vastly expressed in obese subjects as well as up-regulated in hypertensive patients, C1QTNF1 is identified molecular link between obesity and hypertension. C1QTNF1 expression may be associated with a low-grade chronic inflammation status in adipose tissues.

To place an order, please [Click HERE](#).