

## ANXA6 Human

**Description:**ANXA6 produced in E.Coli is a single, non-glycosylated polypeptide chain containing 693 amino acids (1-673 a.a.) and having a molecular mass of 78 kDa. ANXA6 is fused to a 20 amino acid His-Tag at N-terminus and purified by proprietary chromatographic techniques.

**Catalog #:**PRPS-774

For research use only.

**Synonyms:**ANX6, ANXA6, Annexin-VI, Lipocortin-6, Lipocortin-VI, Annexin-6, CPB-II, Annexin A6, p68, p70, Protein III, Chromobindin-20, 67 kDa calelectrin, Calphobindin-II, CBP68.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile filtered colorless solution.

**Amino Acid Sequence:**MGSSHHHHHH SSGLVPRGSH MAKPAQGAKY RGSIHDFPGF  
DPNQDAEALY TAMKGFSDK EAILDIITSR SNRQRQEVQC SYKSLYGKDL IADLYELTG  
KFERLIVGLM RPPAYCDAKE IKDAISGIGT DEKCLIEILA SRTNEQMHQL VAAYKDAYER  
DLEADIIGDT SGHFQKMLVV LLQGTREEDD VVSEDLVQQD VQDLYEAGEL KWGTDEAQFI  
YILGNRSKQH LR

**Purity:**Greater than 90% as determined by SDS-PAGE.

**Formulation:**

The ANXA6 protein solution contains 20mM Tris-HCl, pH-8, 0.1M NaCl and 10% glycerol.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

ANXA6 is part of the family of calcium-dependent membrane and phospholipid binding proteins. ANXA6 mediates the endosome aggregation and vesicle fusion in secreting epithelia during exocytosis. ANXA6 correlates with CD21. ANXA6 regulates the release of Ca(2+) from intracellular stores. ANXA6 is differentially expressed in the lumbar spinal cord from rats submitted to peripheral lesion during neonatal period.

**To place an order, please [Click HERE](#).**