

## PEX19 Human

**Description:** PEX19 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 316 amino acids (1-296) and having a molecular mass of 34.6 kDa. The PEX19 is fused to a 20 amino acid His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: PRPS-932

For research use only.

**Synonyms:** Peroxisomal biogenesis factor 19, Peroxisomal farnesylated protein, HK33, 33kDa housekeeping protein, PXF, PMP1, PXMP1, Peroxin-19, D1S2223E, FLJ55296.

**Source:** E.coli.

**Physical Appearance:** Sterile Filtered colorless solution.

**Amino Acid Sequence:** MGSSHHHHHH SSGLVPRGSH MAAEEGCSV GAEADRELEE  
LLESALDDFD KAKPSPAPPS TTTAPDASGP QKRSPGDTAK DALFASQEKF FQELFDSELA  
SQATAFEKA MKELAEIEPH LVEQFQKLSE AAGRVGSDMT SQQEFTSCLK ETLGLAKNA  
TDLQNSSMSE EELTKAMEGL GMDEGDGEGN ILPIMQSIMQ NLLSKDVLYP SLKEITEKYP  
EWLQSHRESL PP

**Purity:** Greater than 90% as determined by SDS-PAGE.

### Formulation:

The PEX19 solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0) and 10% glycerol.

### Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

### Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Introduction:

PEX19 is vital for peroxisome biogenesis in *Saccharomyces cerevisiae*. PEX19 is limited to the outer surface of peroxisomes in liver cells and is has a part in the early stage of peroxisome membrane assembly, prior to the import of matrix protein.

**To place an order, please [Click HERE](#).**