

Rubella Capsid

Description:The E.Coli derived recombinant protein contains the Rubella Virus Capsid C regions, 1-123 amino acids and fused to a GST-Tag at C-terminus.

Catalog #:RUPS-300

Purity:Rubella protein is >95% pure as determined by 10% PAGE (coomassie staining).

For research use only.

Purification Method:

Rubella protein was purified by proprietary chromatographic technique.

Specificity:

Immunoreactive with sera of Rubella Virus infected individuals.

Formulation:

20mM imidazol, 8M urea and 0.3M NaCl.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Rubella virus is an enveloped positive-strand RNA virus of the family TOGAVIRIDAE: Virions are composed of three structural proteins: a capsid and two membrane-spanning glycoproteins, E2 and E1. During virus assembly, the capsid interacts with genomic RNA to form nucleocapsids.

Storage:

Rubella Capsid protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

To place an order, please [Click HERE](#).